




DEGASi®

Product Program
Degassing & Debubbling

NO TROUBLES WITH BUBBLES®



**BIOTECH
FLUIDICS**



BIOTECH FLUIDICS

Degasser selection guide

FLOW CAPACITY TO AVOID FORMATION OF BUBBLES

Bubbles are detrimental to any fluidic system. Negative pressure transients, heating or mixing solvents reduces gas solubility and causes off-gassing which can lead to bubble formation. Removing up to 50% of dissolved gases with a Biotech DEGASi® in-line degasser effectively prevents bubble formation even in your most demanding applications. Critical factors for choosing the right degasser include flow rate and composition of the liquid. Optimal flow rates to achieve 50% removal are illustrated below. The Constant Vacuum Degassing Control System makes the DEGASi® virtually silent.


BIOCOMPATIBLE FLOW PATH

METAL-FREE

Suitable for proteins, DNA and other biological materials that may adhere to metal surfaces. Not to be used with alkanes like hexane.



Suitable for proteins, DNA and other biological materials that may adhere to metal surfaces. Not to be used with alkanes like hexane.




Flow rate in mL/min

HIGH FLOW RATES FOR WATER-BASED SOLUTIONS
DEGASi® High Flow 1000 with ultra-high degassing efficiency works at flow rates up to 1000 ml/min. Designed for water and water-based solutions like buffers. Can be used with up to 5% water soluble organic solvents.

GPC-COMPLIANT FLOW PATH FOR ANY ORGANIC SOLVENT
GPC-compatible chambers tolerate all common organic solvents including alkanes (hexane, heptane, pentane) and are the recommended choice for DCM, THF, acetone etc. Applications involving Hexafluoroisopropanol (HFIP) require DEGASi® HFIP.

Product details on the back

Let us help you find the right degasser.
biotechfluidics.com



GETTING RID OF BUBBLES - HOW IT WORKS

Biotech DEGASi® degassers use a Systec® AF gas permeable membrane with a Constant Vacuum Degassing Control System to extract dissolved gasses out of fluid. The extracted gasses are evacuated through an exhaust outlet.

DEGASi® product	50% de-gassing mL/min	Max. flow rate mL/min	ID mm	Swept volume ul	Max. fluid pressure bar (psi)	no. of channels	Solvent compatibility			Part number x= no of channels
							Bio-comp.	Alkanes	HFIP	
DEGASi® Plus series										
Nano	0.1	0.3	0.30	25	1.8 (27)	1-6	Yes	No	No	0003-635X-N
Micro	0.5	10	1.14	100	5 (70)	1-6	Yes	No	No	0003-635X-S
Classic	3	10	1.14	480	5 (70)	1-6	Yes	No	No	0003-635X-A
GPC	3	10	1.14	480	5 (70)	1-6	No	Yes	No	0003-662X-B
Semi-Prep	6	10	1.14	925	5 (70)	1-6	Yes	No	No	0003-635X-L
GPC Semi-Prep	6	10	1.14	925	5 (70)	1-6	No	Yes	No	0003-662X-L
DEGASi® Integration series										
Nano	0.1	0.3	0.30	25	1.8 (27)	1-5	Yes	No	No	0001-635X-BN
Micro	0.5	10	1.14	100	5 (70)	1-5	Yes	No	No	0001-635X-BS
Classic	3	10	1.14	480	5 (70)	1-5	Yes	No	No	0001-635X-BA
GPC	3	10	1.14	480	5 (70)	1-5	No	Yes	No	0001-662X-B
Semi-Prep	6	10	1.14	925	5 (70)	1-5	Yes	No	No	0001-635X-BL
GPC Semi-Prep	6	10	1.14	925	5 (70)	1-5	No	Yes	No	0001-662X-BL
HFIP	3	10	1.14	480	5 (70)	1-5	No	Yes	Yes	0001-662X-HFIP
HFIP Semi-Prep	6	10	1.14	925	5 (70)	1-5	No	Yes	Yes	0001-662X-HFIP-L
DEGASi® BioCompact										
BioCompact	2	4	0.89	285	3 (50)	1-6	Yes	Yes	No	0004-X285
DEGASi® Preparative series										
Prep	15	50	1.92	5.4 mL	5 (70)	2 or 4	No	Yes	No	0001-X053
Prep+	75	100	NA	22 mL	3 (50)	1, 2 or 4	No	Yes	No	0001-0X20
Prep+ Large*/XL	100	200	NA	22 mL	3 (50)	1, 2 or 4	No	Yes	No	0001-0X20-L or XL
Prep+ AUX**	75/100	100/200	NA	22 mL	3 (50)	1, 2 or 4	No	Yes	No	0001-0X20-SH
High Flow 1000***	1000	2000	NA	200 mL	1 (14)	1	Yes	No	No	0003-HF.1000

DEGASi® plus



DEGASi® Integration



DEGASi® HFIP



DEGASi® BioCompact



DEGASi® Prep



Fluidic connectors are 1/4"-28 female flatbottom unless otherwise stated.
DEGASi® Nano; 10-32 coned port.

*Large; 5/16"-24 flat bottom female connectors for 3/16" tubing.
XLarge; 1/2"-20 flat bottom female connectors for 1/4" tubing.

**To be complimented with any extra degassing chamber from DEGASi Plus (Micro, Classic, GPC, Semi-Prep or GPC Semi-Prep).

*** Silicone degassing membrane only for water-based solutions and buffers. Barbed end connectors for ID 8 mm tubing.

For ordering and technical support, please contact:

Europe: Biotech Fluidics AB | Tel:+46 (0)300 56 91 80 | info@biotechfluidics.com

USA: Biotech USA LLC | Tel: +1 612 703 5718 | sales@biotechfluidics.com

Japan: BioNik Inc. | Tel: +81 545 38 9125 | info@bionikinc.com

www.biotechfluidics.com



**BIOTECH
FLUIDICS**